

Easy disinfection Washable under running water without detergents





DRAMIŃSKI S.A. Wiktora Steffena 21, Sząbruk 11-036 Gietrzwałd, Poland

↓ +48 89 675 26 00 @ agri@draminski.com

www.draminski.com







MASTITIS DETECTORS

ESSENTIAL EQUIPMENT FOR THE DAILY WORK OF COW, SHEEP AND GOAT FARMERS.

www.draminski.com



OK ESC A



MASTITIS DETECTORS

Essential equipment for the daily work of cow, goat and sheep farmers. Ergonomic models used in forestripping milk from each quarter allow for early diagnosis of udder health.



Ideal equipment for monitoring udder health during the hygienic milking process.

- Enables detection of sub-clinical, asymptomatic inflammation of the udder
- Supports the everyday prevention of mastitis
- Prevents culling cows due to mastitis
- Provides knowledge about the health of breeding animals
- Improves production efficiency and helps to maintain milk quality






MD4X4Q2

MD4Q2

INTERPRETATION OF RESULTS



Left-back quarter: Subclinical inflammation; there is a high risk of disease development and serious inflammation. As a result, medical treatment must be undertaken.

Basic Mode – prompt data interpretation, provides on-site measurement and diagnostics.

Advanced Mode – has all the functions of the basic mode and activates additional device functions, including the option to save readings in the memory and transfer them to a computer for later interpretation.

MAIN FEATURES:

- Waterproof and easy to disinfect
- Intuitive menu
- Large LED display
- Memory for 250 animals and 200,000 measurements with date and time
- Battery powered
- Two-year guarantee

ADDITIONAL FUNCTIONS:

- Real-time clock
- Display mode selection
- Measurement setting (2 or 4 teats)
- Automatic shut-down
- Language selection



J



Operating mode selection Choose function from the menu



PC connection Software update and data transfer